

## **Vibrio spp. Immuno-fluorescence kit**

art. #: ID12322

size: kit for 100/200/400 assays

### **General information *Vibrio spp.***

Several *Vibrio* species are pathogenic to humans. Most disease-causing strains are associated with gastroenteritis, but can also infect open wounds and cause sepsis. It can be carried by numerous sea-living animals, such as crabs, oysters or prawns, and has been known to cause fatal infections in immune-compromised humans during exposure. Examples of pathogenic *Vibrio* species include *V. cholerae*, *V. parahaemolyticus*, and *V. vulnificus*. *Vibrio cholerae* is the causative agent of Cholera infections. Cholera infections are most commonly acquired from drinking water in which *V. cholerae* is found naturally or into which it has been introduced by fecal contamination. *V. cholerae* colonizes the small intestine and releases enterotoxins. As a consequence of improved sewage and water treatment, there is a low occurrence of *V. cholerae* in the Western civilizations, but *V. cholera* still possess a serious threat in third world counties and after flooding of civilized areas. Over 90% of the Cholera cases occurring in the Western civilizations are the result of travel to a country where this bacterium is prevalent. *V. vulnificus* and *V. parahaemolyticus* are responsible for the vast majority of seafood-related infections. In healthy people, ingestion of *V. vulnificus* and *V. parahaemolyticus* can cause vomiting, diarrhea, abdominal pain and skin infections. In immune-compromised patients, especially those with chronic liver disease, are at high risk for *V. vulnificus* when they eat raw seafood, particularly oysters. A recent study showed that people with these pre-existing medical conditions were 80 times more likely to develop *V. vulnificus* and *V. parahaemolyticus* bloodstream infections than were healthy people. *V. vulnificus* and *V. Parahaemolyticus* are frequently isolated from oysters and other shellfish in warm coastal waters during the summer months. Since it is naturally found in warm marine waters, people with open wounds can be exposed to *V. vulnificus* and *V. parahaemolyticus* through direct contact with seawater.

### **Applications**

Excellent suitable for immunofluorescence staining procedures on glass, other applications or platforms are not tested but should not be excluded.

### **Contents:**

<b>Content</b>	<b>Format</b>	<b>Use</b>	<b>Store at</b>
Wash Buffer 1	Liquid	Ready to use	Room temperature
Blocking Buffer	Liquid	Ready to use	4°C, before use warm up to room temperature
Wash Buffer 2	Liquid	Ready to use	Room temperature

### **Vibrio spp. antibody**

**Clonality** : polyclonal  
**Immunogen** : whole cells of various strains of *Vibrio*  
**Host animal** : rabbit  
**Conjugation** : fluorescein analogue  
**Purification** : IgG purified  
**Format** : lyophilized

---

**Stabilizer and preservative**

IgG free bovine serum albumin (BSA) is added as a protein stabilizer. 0,02% sodium azide is added as a preservative. Non-sterile.

**Antibody concentration**

The concentration of affinity purified antibody is 0,5 mg as determined by UV absorbance at 280nm. Upon rehydration with water, the solution will contain 1% BSA, 100mM phosphate, 150mM sodium chloride, 0,02% sodium azide, pH 7,4. The conjugate will be at a concentration of 0,5 mg/ml.

**Rehydration**

Rehydrate with 1 ml reagent quality water. Rotate the vial until the lyophilized pellet is totally dissolved.

**Use**

Dilute to the desired concentration with blocking buffer immediately before use, mix thoroughly. This working solution is not recommended for long term storage.

**Storage**

Store at 2-8°C until rehydration, rehydrated antibody may be stored for up to one week at 2-8°C, thereafter it should be stored at -20°C. Avoid multiple freeze thaw steps. When aliquoting, store product in volumes greater than 50µl. Variations in temperature due to freeze cycles may cause loss of activity when rehydrated product is stored frozen in aliquots less than 50µl.

**Specificity**

This antibody broadly reacts to *Vibrio* species, including *Vibrio cholera O1*, *Vibrio cholera non O1*, *Vibrio parahaemolyticus* and *Vibrio vulnificus*.

**Excitation/emission values**

Fluorescein analogue is excited at 493 nm (in PBS) and emits at 518 nm (in PBS).

**Contact information**

If you have any questions about this product, please contact us at [Sales@innosieve.com](mailto:Sales@innosieve.com) or call us at (+31)-646717500.

### Protocol in eppendorf tube:

#### Notes before starting:

- All centrifugation steps are performed for 2 minutes at 14.000 RCF (relative centrifugal force)
- Not all the supernatant is removed to prevent loss of the pellet

#### Method:

1. Pipette 500µl sample in an eppendorf tube  
*Remark: if the amount of sample is less, add Wash Buffer 1 to a final volume of 500µl*
2. Vortex the sample and centrifuge
3. Remove 450µl of the supernatant without disturbing the pellet
4. Resuspend the pellet in the remaining supernatant
5. Add 500µl Blocking Buffer, vortex and centrifuge, remove 500µl of the supernatant
6. Resuspend the pellet in the remaining supernatant
7. Add 500µl Wash Buffer 1, vortex and centrifuge, remove 500µl of the supernatant
8. Resuspend the pellet in the remaining supernatant
9. Prepare the antibody solution, dilute for one sample 10µl antibody stock in 90µl Blocking Buffer, mix by pipetting up and down
10. Add 100µl prepared antibody working solution, vortex and incubate at room temperature for 10 minutes
11. Add 400µl Wash Buffer 1, vortex and centrifuge, remove 500µl of the supernatant
12. Resuspend the pellet in the remaining supernatant
13. Add 500µl Wash Buffer 1, vortex and centrifuge, remove 500µl of the supernatant
14. Resuspend the pellet in the remaining supernatant  
*Optional: in case higher stringency is required add 500µl Wash Buffer 2, vortex and centrifuge. Discard the supernatant and resuspend the pellet in the remaining supernatant. This step can be repeated. When applying this stringency step, finish by applying Wash Buffer 1 to resuspend the cells.*
15. Add 5µl of the bacterial suspension onto a glass slide and allow to air dry in the dark
16. Heat fix the sample by passing the glass slide the flame for 3 or 4 times
17. Add Mounting Medium and cover glass
18. Analyze the sample

## Protocol on glass slide

### Notes before starting:

- Any type of sample can be used, preferably suspended cells in buffer

### Method 1:

1. Add the sample onto a glass slide and allow to air dry in the dark  
*Note: in most cases a sample volume of 5 - 50  $\mu$ l is recommended*
2. Heat fix the sample by passing the glass slide the flame for 3 or 4 times
3. Add 50 $\mu$ l blocking buffer, allow to stand 1 minute, rinse gently with 500 $\mu$ l wash buffer 1
4. Prepare the antibody solution, dilute for one sample 5 $\mu$ l antibody stock in 45 $\mu$ l blocking buffer, mix by pipetting up and down
5. Add 50 $\mu$ l diluted antibody, pipette gently up and down 8 times and incubate at RT for 10 minutes  
*Note: depending on the sample the incubation time can be elongated to improve the signal*
6. Rinse twice gently with 500 $\mu$ l wash buffer 1  
*Optional: in case higher stringency is required rinse gently with 500 $\mu$ l Wash Buffer 2. This step can be repeated. When applying this stringency step, finish with one wash step using Wash Buffer 1*
7. Allow to air dry the sample in the dark
8. Add Mounting Medium and cover glass
9. Analyze the sample

## Protocol on glass slide

### Notes before starting:

- Any type of sample can be used, preferably suspended cells in buffer

### Method 2:

1. Add the sample onto a glass slide and allow to air dry in the dark  
*Note: in most cases a sample volume of 5 - 50  $\mu$ l is recommended*
2. Heat fix the sample by passing the glass slide the flame for 3 or 4 times
3. Rinse twice gently with 500 $\mu$ l MilliQ
4. Add 50 $\mu$ l Wash Buffer 3, allow to stand 2 minutes, rinse gently with 500 $\mu$ l Wash Buffer 3
5. Rinse gently with 500 $\mu$ l blocking buffer
6. Prepare the antibody solution, dilute for one sample 5 $\mu$ l antibody stock in 45 $\mu$ l Blocking Buffer, mix by pipetting up and down
7. Add 50 $\mu$ l diluted antibody, pipette gently up and down 8 times and incubate at RT for 10 minutes
8. *Note: depending on the sample the incubation time can be elongated to improve the signal*
9. Rinse thrice gently with 500 $\mu$ l Wash Buffer 1  
*Optional: in case higher stringency is required rinse gently with 500 $\mu$ l Wash Buffer 2. This step can be repeated. When applying this stringency step, finish with one wash step using Wash Buffer 1*
10. Allow to air dry the sample in the dark
11. Add Mounting Medium and cover glass
12. Analyze the sample